



# GenoSensor

## EduPrimer™ DNA Profiling Kit

### Brief Protocol

**EduPrimer™ DNA Profiling Kit is specifically designed for exposing novice students to PCR principles and technique. The kit is simple to use. The whole lab can be done within 3 hours.**

#### DNA Preparation ~ 7 min

1. Add 200µL of **Solution A** to a 1.5mL microcentrifuge tube.
2. Thoroughly swab inside cheek with provided swab and put it into **Solution A**.
3. Vortex the sample for 10 seconds, then heat in pre-heated 95°C heat block for 5 minutes
4. Briefly spin sample down in microcentrifuge. Remove the swab with tweezers.
5. Add 20µl **Solution B** to the sample tube. Vortex or invert to mix for at least 10 seconds.
6. Spin sample for 1 minute at 12,000rpm
7. Use 10 µl of supernatant as DNA template for PCR.

#### PCR Reaction Mixture ~ 5 min

Mix the following reagents into a standard PCR tube:

2X PCR Master Mix	10µl
<u>Genomic DNA Template</u>	<u>10µl</u>
Volume total = 20µl	

#### PCR Parameters ~ 78 min

1. 94°C – 2 minutes
2. 94°C denaturing – 20 seconds}
3. 58°C annealing – 20 seconds} **repeat steps 2, 3, & 4 for 40 cycles**
4. 72°C extension – 20 seconds}
5. 72°C – 5 minutes
6. 4°C – finished / hold



#### Agarose Gel Electrophoresis ~30 min

- Pour 1% agarose gel, using your preferred staining method\*.
- Use 10 µL of PCR product in each well to visualize on gel.
- Run at ~100V for 10-20 minutes and stop before loading dye has run off gel
- Visualize and record the results manually or by photography
  - Larger expected band – **400bp** (Alu element inserted)
  - Smaller expected band – **100bp** (no Alu insert)

\*It's recommended that an in-gel stain visualized by UV light is used as they provide the sensitivity necessary for the smaller fragments amplified in this experiment.

#### Additional Required Materials

Thermal cycler, Heat block, Microcentrifuge, Micropipettes, Pipette tips, PCR tubes, Gel electrophoresis apparatus

(Full protocol for students also available on [our website](#))